

当显示屏在每个校准点显示 pH 7 时该怎么办？

回复：

须先跟您确认产品序号来确认发货时间。若该物品已经发货一年多，我们假设原因是因为电极老化。

尽管如此，您仍然可以尝试通过将电极浸入水龙头 10-15 分钟来清洁探头，然后用蒸馏水彻底冲洗，看看它是否仍然可用的。

Question:

What to do when display shows pH 7 at every calibration points?

Reply:

If the item already had been shipped out for more than one year, we assume that the reason is because of electrode aging.

Nevertheless, you can still try to clean the probe by immersing the electrode in tap water for 10-15 minutes, then rinse thoroughly with distilled water, to see if it's still available.